

ROTUNDA HOSPITAL DEPARTMENT OF LABORATORY MEDICINE
POCT Active Test Repertoire Table LF-GEN-0076 Ed 02

Please note: Printed reference ranges from the blood gas analysers, and ranges displayed on the screen of the Hb and glucose meters are not configurable to pregnancy/non pregnancy status. Samples processed in **syringe** mode on the blood gas analysers will print **adult non pregnant** reference ranges. Samples processed in **capillary mode** on the blood gas analysers will print **neonatal reference ranges**. Correct sample type should be selected e.g. arterial, venous, cord, fetal scalp should be selected to ensure appropriate ranges apply. Ranges displayed on the glucose and Hb meters are for adult non pregnant patients. See table below for a more comprehensive list of ranges related to age, pregnancy status and sample type

Results should be interpreted in relation to agreed clinical decision limits. Where a POC result warrants clinical intervention, or where the care giver is not reassured by the result, a sample should be sent to the laboratory for confirmation

TEST	SAMPLE TYPE	SPECIAL PRECAUTIONS	REFERENCE RANGE	Reference	INAB ACCREDITED TEST REG NO 208MT
pH	Venous, Arterial, Capillary, Cord vein/artery, fetal scalp Via Syringe, capillary tube, as appropriate	Ensure no air pockets	Adult Whole Blood Arterial 7.35-7.45 Venous 7.32-7.43	1	Yes
		Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.	Cord Blood Arterial 7.18-7.38 Venous 7.25-7.45	1	
		Capillary tube samples- Ensure no air pockets Fetal Scalp tube-Ensure the end of tube that has touched baby's scalp is not presented to analyser	Fetal Scalp >= 7.25 Normal 7.21-7.24 Borderline <= 7.2 Abnormal Result	2	
pCO2		Ensure no air pockets Syringe Sample-mix well by inverting several times and	Adult Whole Blood, Arterial 4.3-6.0 kPa	1	Yes

	<p>Venous, Arterial, Capillary, Cord vein or artery, Fetal scalp.</p> <p>Via Syringe or capillary tube, as appropriate</p>	<p>rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe</p>	<p>Infant WB Arterial 3.6-5.5kPa</p>	1	
pO2	<p>Venous, Arterial, Capillary, Cord vein or artery, Fetal scalp.</p> <p>Via Syringe , capillary tube, as appropriate</p>	<p>Ensure no air pockets</p> <p>Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.</p> <p>Capillary tube samples- Ensure no air pockets</p> <p>Fetal Scalp -Ensure the end of tube that has touched baby's scalp is not presented to analyser.</p>	<p>2 days-60 years Whole Blood arterial 11.04-14.4 kPa</p> <p>Cord Blood Arterial 0.8-4.1kPa Venous 2.3-5.5kPa</p>	1 1	Yes
Sodium	<p>Venous, Arterial, Capillary. Via syringe or capillary tube, as appropriate.</p>	<p>Ensure no air pockets in sample</p> <p>Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.</p>	<p>Pregnant, Li-Hep Plasma 133-143 mmol/L</p> <p>Non pregnant, Li-Hep Plasma 133-146 mmol/L</p> <p>Infants, Li-Hep Plasma 133-144 mmol/L</p>	3 4 5	Yes

Potassium	Venous, Arterial, Capillary Via Syringe or capillary tube, as appropriate	Ensure no air pockets Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling. Vacutainer-mix by inverting and rolling several times. Sample via sample adaptor	Pregnant, Li-Hep Plasma 3.5-5.3mmol/L	3	Yes
			Non Pregnant, Li-Hep Plasma 3.5-5.3mmol/L	1	
Chloride	Venous, Arterial, Capillary, Via Syringe or capillary tube as appropriate	Ensure no air pockets Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.	Neonate, Li-hep Plasma 3.5-5.5mmol/L	5	Yes
			Pregnant, Li-hep Plasma 95-108mmol/L	3	
			Non-Pregnant 95-108mmol/L	1	
Neonates, Li-Hep plasma 98-110mmol/L	5				
Glucose-Blood Gas Analyser:	Venous, Arterial, Capillary. Via Syringe or capillary tube,as appropriate	Ensure no air pockets Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.	Pregnant, non pregnant Fasting, neonates Plasma 3.0-5.1mmol/L	6	Yes

Glucose Glucose Meter	Wholeblood- Fingerprick/Heel Prick (results are plasma equivalent)	Do not wipe puncture site with alcohol wipe before sampling Use sterile water to clean Sample site	Pregnant, non pregnant fasting, & neonates Plasma 3.0-5.1mmol/L	6	Yes
Lactate Lactate Contd..	Venous, Arterial, Capillary, Via Syringe or capillary tube , as appropriate	Ensure no air pockets Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling. Capillary tube samples-Ensure No air pockets	Adult WB heparin Venous 0.6-1.4mmol/L Arterial 0.4-0.8mmol/L At Birth Umbilical Artery 1.6-5.5 mmol/L Umbilical Vein 1.2-5.0mmol/L Neonates <2.9mmol/L	1 7 8	Yes Yes
Ionised Calcium	Venous, Arterial, Capillary Via Syringe or capillary tube as appropriate	Ensure no air pockets Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling. Capillary tube samples-Ensure no air pockets	Adult. WB (Heparin) 1.18 -1.32mmol/L Infants 0-1month 1.0-1.5mmol/L Infants 1-6 months 0.95-1.5mmol/L	9	Yes
Bilirubin	Venous, Arterial, Capillary	Ensure no air pockets Syringe Sample-mix well by			

Hb-Blood Gas Analyser	Via Syringe or capillary tube as appropriate	<p>inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.</p> <p>Capillary tube samples- Ensure no air pockets</p>	<p>Non pregnant, Li-hep Plasma <21umol/L</p> <p>Pregnant, Li-hep Plasma <21umol/L</p> <p>Neonates, Li-hep Plasma <17umol/L</p>	10 10 10	Yes
	Venous, Arterial, Capillary, Via Syringe or Capillary tube, as appropriate	<p>Ensure no air pockets</p> <p>Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.</p>	<p>Pregnant, Whole Blood EDTA 100-160g/L</p> <p>Non-Pregnant, whole blood EDTA 120-160g/L</p> <p>Infant Whole Blood EDTA</p> <p>Day 0-2 : 135-195g/L Day 2-4 : 145-255g/L Day 5-8: 135-215g/L Day 9-21: 125-205g/L Day 22-35:110-180g/L Week 6-9 : 90-140g/L Week 10-Month 18: 110-135g/L</p>	11 11 11 11	
Hb-Hemocue Hb meter	Finger Prick/heel Prick	Expel first 2-3 drops from puncture site before sampling.	Pregnant, EDTA Whole Blood 100-160g/L	11	Yes

		Wipe cuvette with tissue before sampling-ensure no blood is drawn out of cuvette. Ensure cuvette is completely Filled and free from air bubbles	Non-Pregnant, EDTA Whole Blood 120-160g/L Infant EDTA Whole Blood Day 0-2 : 135-195g/L Day 2-4 : 145-255g/L Day 5-8: 135-215g/L Day 9-21: 125-205g/L Day 22-35:110-180g/L Week 6-9 : 90-140g/L Week 10-Month 18: 110-135g/L	11 11	
MethHb	Venous, Arterial. Capillary. Via Syringe or capillary tube as appropriate	Ensure no air pockets Syringe Sample-mix well by inverting several times and rolling between palms. Expel first drop of blood onto tissue. Ensure blood is at top of syringe before sampling.	Pregnant, Non Pregnant, Infant, Whole Blood *0.0-1.5%	1	Yes
Covid POC Abbott ID Now	Nasal Swab	Inset swab until you meet resistance. Touch swab off nasal wall and rotate several times, repeat with same swab on both nostrils	N/A		No
Pregnancy testing-urine	Urine	First morning urine recommended as is contains highest conc. of hcg. Not less than 1ml of sample.	N/A Sensitivity of 25mIU/ml	12	No

References

1. Teitz Textbook of Clinical Chemistry and Molecular Diagnostics 5th Ed. 2012
2. Rotunda Hospital Guidelines on Intrapartum Fetal Surveillance PPGC-OBC-21 ; Nice 2007
3. In House data, Rotunda Hospital. Biochemistry central lab reference ranges REF-LP-GEN-0007
4. Pathology Harmonisation UK2011. Biochemistry central lab reference ranges REF-LP-GEN-0007
5. In house data from Temple St Children's Hospital and OLHSC, Crumlin. Biochemistry central lab reference ranges REF-LP-GEN-0007
6. HSE Guidelines 2007. Biochemistry central lab reference ranges REF-LP-GEN-0007.
7. Olesen H. Kompendium 2000 - Kompendium i laboratoriemedicin. <http://www.dskb.dk/default.asp?id=327> as accessed 21-02-2014.
8. Heil W, Koberstein R, Zawata B. Reference ranges for adults and children; preanalytical considerations. Wuppertal: Roche Diagnostics GmbH, 2001.
9. Soldin SJ, Wong EC, Brugnara C et al. Pediatric reference intervals. 7th ed. Washington DC: AACC Press, 2011.
10. Biochemistry central lab reference ranges REF-LP-GEN-0007
11. Sourced Internally-refer to Haematology Laboratory Reference ranges for source
12. SureScreen Diagnostics Kit Insert